ENCISION recommends placing this document in the Instructions for Use/Care Section of your AEM Monitor Operator/Service Manual.

Device Description

AEM Universal Adapter

The **ES9015** Universal Adapter is intended to be used only with Encision's AEM Monitoring System. The Universal Adapter provides interface between the AEM Handswitching Handpiece (ES1300 series), the AEM Cord (ES4007 or ES4107) for footswitching instruments, the AEM Monitor, and compatible electro-surgical generators.


NOTE

- This product is rated to 9000 Vp-p. Limit electro-surgical generator power setting to 80 Watts or lower (60 Watts for the Commed Aspen Excalibur spray mode). Higher settings may result in spurious insulation failure alarms and/or insulation breakdown. See Encision AEM Monitor Operator/Service Manual for list of compatible electro-surgical generators.
- This product is not intended to be located in the sterile field. It is not supplied sterile. Steam sterilization will damage the device.

Use with Monopolar Electro-surgery

AEM instruments, in conjunction with an AEM Monitor properly connected to the electro-surgical generator (ESU), continuously monitor and dynamically manage "stray energy" (insulation failure and capacitive coupling) in zones 2 & 3, which are likely out of the surgeon's field of view.

AEM shielding does not cover zone 1, which the surgeon should keep in view during instrument activation. As in all applications, "misapplied" electro-surgical energy remains the responsibility of the attending surgeon.

(Note: Zone 1 equals Tip to Shield)
(Note: Zone 3 equals area of Trocar Cannula)

Electro-surgical procedures should be performed only by surgeons having adequate training and familiarity with these techniques and who are also knowledgeable about anatomy and pathology as well as the complications, hazards, risks and benefits of the procedure.

Indications/Intended Use

This accessory is for AEM instruments which are electro-surgical cutting and coagulation accessories intended to remove tissue and control bleeding by use of high-frequency electrical current.

AEM instruments incorporate the use of AEM technology and are intended for use with the AEM Monitoring System and electro-surgical generators having compatibility with the AEM Monitor.

Contraindications

These instruments are not intended for use when electro-surgical techniques are contraindicated.

Instructions For Use

Prior to Use

Thoroughly read these instructions and the instructions in the AEM Monitor Operator/Service Manual.

Inspect the instrument for proper assembly and function.

AEM System Setup

See laminated Setup Sheet (00701) when using the ES9005 series AEM Cord Adapter and (02678) when using the ES9015 Universal Adapter.

See laminated Troubleshooting Guide (02678) for additional information.

WARNING

- Pulling or tugging cords can result in breaking of internal conductors, causing intermittent alarms or sparking and burning of insulation during use. The AEM Monitor is not designed to detect or prevent an arc in the event of a broken active wire.
- Keep electrical connections dry while in use to prevent potential conduction of HF current to the user.
- Damaged external insulation AND incorrect setup of the AEM Monitor may result in a risk of unintended patient burn, shock or fire hazard. Do not use product having damaged insulation.

CAUTION

- Good operating room practice suggests that connections of accessories to electro-surgical generators be made only while the generator is Off or on Standby.
- Use these instruments only in conditions that assure adequate visualization to minimize risk of misapplied electro-surgical energy.
- Do not activate the ESU if the monitor's green Ready light is off. No Ready light means there is a setup error or equipment malfunction.
- Keep ESU power setting as low as possible for the intended purpose to minimize unintended burns.
- Damaged internal insulation of the cord and/or instrument, or loss of shield continuity, may cause ESU return pad alarms triggered by the AEM Monitor's Fault Indicators. For maximum patient safety, discontinue use of the instrument if this occurs.
- A singular AEM instrument must be the sole conductor of energy to tissue. Do not conduct energy by touching an AEM instrument to a second instrument contacting tissue. The second device will not be protected from capacitive coupling and insulation failure.
- Keep electro-surgical instruments away from the patient and operative field when not in use. Accidental activation can result in unintended injury to the patient.
- See electro-surgical generator manual and AEM Monitor Operator/Service Manual for precautions concerning the general application of electro-surgical equipment.

Universal Adapter Installation

A Philips Head screw driver is required to install the Universal adaptor on the monitor.

Each Universal Adapter kit contains:


- (1) Universal Adapter
- (1) 8-32 x 3/4" long Philips head screw
- (1) 10-24 x 1/2" long Philips head screw
- (1) Troubleshooting Guide (02678)
- (1) 6-32 x 1/2" long Phillips head screw
- (1) #6 washer
- (1) #8 washer

Optional – Remove the Inhibit Adapter Extension:


- (1) 4-40 x 3/4" long hex screw
- (1) Allen Wrench

Use the allen wrench to unscrew the existing 4-40 x 1" screw. Discard inhibit adapter extension. Use the allen wrench to attach the 4-40 x 3/4" screw to the inhibit adapter.

1. Remove and discard the front 8-32 UNC screw from right side of the AEM Monitor enclosure and the left-hand 6-32 UNC screw from the AEM Cord receptacle on the front panel of the AEM Monitor.

2. Plug the Inhibit Adapter into the BLUE receptacle on the Universal Adapter.
3. Push the Inhibit Adapter cord through the vertical groove on the right side of the Universal Adapter.
4. Secure the Inhibit Adapter to the Universal Adapter using the 10-24 UNC screw provided.
5. Plug the Universal Adapter into the AEM Cord receptacle of the AEM Monitor

6. Firmly attach the Universal Adapter to the front of the AEM Monitor using the 6-32 screw of the AEM and washer provided. Attach it to the side of the Monitor using the 8-32 screw and washer provided.
7. Attach the Troubleshooting Guide to the key ring on the side of the Universal Adapter.

CAUTION


- Do not use the AEM Monitor without both the GREEN Inhibit Adapter plugged into the Universal Adapter and the Universal Adapter's BLUE-banded Return Electrode cord plug fully plugged into the ESU socket. This is the means by which the AEM Monitor signals the ESU to stop output when insulation or connection failure is detected.

AEM Monitor setup

1. Turn on AEM Monitor. Verify all visual indicators briefly illuminate, and audible beeps are heard.
2. Verify the Universal Adapter is properly installed on the AEM Monitor.
3. Plug the GREEN-banded cord plug into the ESU Foot Control Receptacle.
4. Plug the YELLOW-banded cord plug into the ESU Hand Control Receptacle.
5. Plug the BLUE-banded cord Plug into the ESU Return Electrode Receptacle.

Instrument setup

1. Plug the Patient Return Electrode Pad into the Inhibit Adapter. Verify that the amber Return Electrode light turns off. The system is now ready for use with bipolar and non-AEM monopolar instruments.
2. For a foot controlled instrument, plug the AEM Cord into the GREEN receptacle on the Universal Adapter, and plug the other end of the cord into the instrument.
3. For a hand controlled instrument, plug the instrument cord into the YELLOW receptacle on the Universal Adapter.
4. After installing all desired instruments, verify that the GREEN Ready light is ON. The AEM Monitoring System is ready for use.

NOTE

- The installation of AEM instruments does not affect the use or performance of non-AEM or Bipolar instruments. The order of installation is not critical.
- Always disconnect cables by pulling on the connector body. Replace adapter if intermittent alarms are experienced during use or setup.

CAUTION

- Do not use the AEM Monitor without both the GREEN Inhibit Adapter plugged into the Universal Adapter and the Universal Adapter's BLUE-banded Return Electrode cord plug fully plugged into the ESU socket. This is the means by which the AEM Monitor signals the ESU to stop output when insulation or connection failure is detected.

Reprocessing

Cleaning

Wipe down the exterior using a mild cleaning solution or disinfectant and a damp cloth to clean the adapter.

Use a cotton swab to clean the internal contacts of any contaminates.

NOTE

- Do not steam sterilize or immerse this product in any liquid. This will damage the device.
- Used instruments are considered medical waste. Dispose of in accordance with local regulations.

End of Life Indicators

Visually inspect and electrically test prior to use. Discontinue use if any of the following are evident:

- Connector: Bent, broken, loose or missing components, or corroded contacts. Damaged active contact (center post) insulation (i.e. cracks, nicks, abrasions, holes, tears, burns or melted areas).
- Cord: Damaged insulation (i.e. Cracks, nicks, abrasions, holes, tears, kinks, bulges, burns or melted areas). Pay special attention to cord/connector interface. Manually flex the cable near the connector ends. Increased flexibility may indicate a conductor wire break, even without obvious insulation damage.
- Electrical: Verify proper operation with the AEM Monitor and electro-surgical generator; the adapter should be removed from service if any malfunction occurs. Refer to "Checking the AEM Monitor" in the AEM Monitor Operator/Service Manual.

Express Warranty

ENCISION hereby warrants to Buyer that products purchased hereunder shall be free from defects in materials and workmanship under normal use and service, as specified in this Instruction for Use/Care, for the period of 6 months from date of purchase. Any evidence of unauthorized repair, modification, or sterilization of this product will void the warranty.

See AEM Monitor Operator/Service Manual for details of Limitations, Disclaimer, and Exclusions.

Return of Used Product

If for any reason this product must be returned to ENCISION, a returned goods authorization is required prior to shipping. Appropriate return instructions may be obtained from ENCISION.

Product

ENCISION reserves the right to amend, modify or to change any product, to introduce new products, to withdraw products and otherwise vary product specifications at any time without notice.

ENCISION® and AEM® are registered trademarks of ENCISION Inc.

Rx ONLY	Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult Instructions for Use
	Latex Free
	Foot Switch
	Hand Switch
	Patient Return Electrode Jack

Made in USA

Manufactured by

ENCISION Inc.
6797 Winchester Circle
Boulder, CO 80301 USA
Ph: 303-444-2600
Fax: 303-444-2693

Authorized Representative

(according to MDD93/42/EEC
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Printed in USA
© Copyright 2010 Encision Inc.

This page intentionally left blank.

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ENCISION® en AEM® zijn gedeponeerde handelsmerken van ENCISION Inc.