# Inhibit Adapter Extension

Instructions For Installation and Use/Care

EN


**Rx ONLY** 

ENCISION recommends placing this document in the Instructions for Use/Care Section of your AEM® Monitor Operator/Service Manual.


## Device Description

### Inhibit Adapter Extension


Inhibit Adapter Extensions are accessories designed for use with the existing Inhibit Adapter of the Encision AEM® Monitor.

**ES9007** Inhibit Adapter Extension  
– for use with Conmed Monitors

**ES9008** Inhibit Adapter Extension  
– for use with monitors needing an extension adapter as recommended by Encision

**ES9009** Inhibit Adapter Extension  
– for use with Olympus Monitors


**NOTE**

- This product is rated to 9000 Vp-p. Limit electro-surgical generator power setting to 80 Watts or lower (60 Watts for the Conmed Aspen Excalibur spray mode). Higher settings may result in spurious insulation failure alarms and/or insulation breakdown. See Encision AEM Monitor Operator/Service Manual for list of compatible electro-surgical generators.
- This product is not intended to be located in the sterile field. It is not supplied sterile. Steam sterilization will damage the device.

### Use with Monopolar Electrosurgery

AEM instruments, in conjunction with an AEM Monitor properly connected to the electro-surgical generator (ESU), continuously monitor and dynamically manage “stray energy” (insulation failure and capacitive coupling) in zones 2 & 3, which are likely out of the surgeon’s field of view.

AEM shielding does not cover zone 1, which the surgeon should keep in view during instrument activation. As in all applications, “misapplied” electro-surgical energy remains the responsibility of the attending surgeon.

(Note: Zone 1 Tip to Shield)  
(Note: Zone 3 equals area of Trocar Cannula)

Laparoscopic procedures should be performed only by surgeons having adequate training and familiarity with laparoscopic techniques and who are also knowledgeable about anatomy and pathology as well as the complications, hazards, risks and benefits of the procedure.

### Indications/Intended Use

These AEM instruments for which this is an accessory incorporate the use of AEM technology and are intended for use in delivering monopolar electro-surgical energy during laparoscopic procedures only.

AEM instruments are intended for use with the AEM Monitoring System and electro-surgical generators having compatibility with the AEM Monitor.

### Contraindications

These instruments are not intended for use when laparoscopic electro-surgical techniques are contraindicated.

## Instructions For Use

### Prior to Use

Thoroughly read these instructions and the instructions in the AEM Monitor Operator/Service Manual.

Inspect the instrument for proper assembly and function (see End of Life Indicators).


### Installation Instructions

All necessary tools and hardware are provided with the Inhibit Adapter Extension. See general illustration below.


1. Use provided Allen wrench to unscrew one of the 3/4 inch hex screws from the inhibit adapter. This screw is no longer needed.

**NOTE**


- Inhibit adapter extension may be removed if the adapter extension is no longer desired. If this is desired, the 1 inch hex screw must be replaced with the original 3/4 inch hex screw.

2. Attach male end of inhibit adapter to female end of inhibit adapter extension.

3. Use Allen wrench to screw in provided 1 inch hex screw. Hand tighten; do not over torque as adapter housing may crack.

4. Perform AEM monitor functions per section 3 of the AEM monitor operator/service manual. Perform normal electrical safety tests as required per your Hospital’s Protocol.

### AEM System Setup

See laminated Setup Sheet (00701) when using the ES9005 series AEM Cord Adapter and (02678) when using the ES9015 Universal Adapter.

#### WARNING

- Pulling or tugging cords can result in breaking of internal conductors, causing intermittent alarms or sparking and burning of insulation during use. The AEM Monitor is not designed to detect or prevent an arc in the event of a broken active wire.
- Damage to connector body and/or insulation may result in shock or fire hazard.
- Keep electrical connections dry while in use to prevent potential conduction of HF current to the user.
- Damaged external insulation AND incorrect setup of the AEM Monitor may result in a risk of unintended patient burn, shock or fire hazard. Do not use product having damaged insulation.

#### CAUTION

- Good operating room practice suggests that connections of accessories to electrosurgical generators be made only while the generator is Off or on Standby.
- Keep ESU power setting as low as possible for the intended purpose to minimize unintended burns.
- A singular AEM instrument must be the sole conductor of energy to tissue. Do not conduct energy by touching an AEM instrument to a second instrument contacting tissue. The second device will not be protected from capacitive coupling and insulation failure.
- Keep electrosurgical instruments away from the patient and operative field when not in use. Accidental activation can result in unintended injury to the patient.
- Always fully connect AEM cord to instrument and adapter. Partial connections may cause arcing, which will damage the product.
- See electrosurgical generator manual and AEM Monitor Operator/Service Manual for precautions concerning the general application of electrosurgical equipment.

### Reprocessing Instructions

#### Manual Cleaning

Wipe down the exterior using a mild cleaning solution or disinfectant and a damp cloth to clean the adapter.

#### Automatic Cleaning

Do not place this product in an automatic washer or immerse this product in any liquid. This will damage the device.

#### Sterilization

Do not sterilize. Exposure to steam sterilization or cold soak sterilization will damage the device.

#### NOTE

- Used instruments are considered medical waste. Dispose of in accordance with local regulations.

### Product Life

The life of this surgical instrument is largely dependent on the care and handling at the point of use or cleaning. For optimal instrument life, protect it from contact / impact with other instruments during decontamination. The number of uses may be reduced by improper handling.

### End of Life Indicators

Visually inspect and electrically test prior to use. Discontinue use if any of the following are evident:

- Connector: Bent, broken, loose or missing components, or corroded contacts. Damaged active contact (center post) insulation (i.e. cracks, nicks, abrasions, holes, tears, burns or melted areas).
- Electrical: Verify proper operation with the Active Electrode Monitoring System and electrosurgical generator; the adapter should be removed from service if any malfunction occurs. Refer to "Checking the Active Electrode Monitoring System" in the AEM Monitor Operator/Service Manual.

### Express Warranty

ENCISION hereby warrants to Buyer that products purchased hereunder shall be free from defects in materials and workmanship under normal use and service, as specified in these Instructions for Use/Care, for the period of one year from date of purchase, or upon reaching an end of life indicator, whichever occurs first. Any evidence of unauthorized modification, repair or sterilization of this device will void the warranty.

See AEM Monitor Operator/Service Manual for details of Limitations, Disclaimer, and Exclusions.

### Return of Used Product

If for any reason this product must be returned to ENCISION, a returned goods authorization is required prior to shipping. Appropriate return instructions may be obtained from ENCISION.

### Product

ENCISION reserves the right to amend, modify or to change any product, to introduce new products, to withdraw products and otherwise vary product specifications at any time without notice.

ENCISION® and AEM® are registered trademarks of ENCISION Inc.

 <b>Rx ONLY</b>	Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult Instructions for Use
	Latex Free

Made in USA

#### **Manufactured by**

ENCISION Inc.  
6797 Winchester Circle  
Boulder, CO 80301 USA

Ph: 303-444-2600  
Fax: 303-444-2693

Applies to ES9007 Only

#### **Authorized Representative**

(according to MDD93/42/EEC  
MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany

Printed in USA  
© Copyright 2010 Encision Inc

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ENCISION® en AEM® zijn gedeponeerde handelsmerken van ENCISION Inc.