# Reusable Articulating Inserts

Instructions For Use/Care

REF ES05XX Series

EN

## Device Description


### Instrument Inserts ES05xx Series

Most inserts fit standard 5.5mm trocar cannulas. (Inserts with a larger or smaller cross-section are identified in the catalog description.)

### Principle of Operation

The front and rear handles open and close the working tip (jaw/ blade insert) of the instrument.

The working tip of the instrument should always be closed when introducing or removing the instrument from the cannula.

### NOTE

- This insert has been validated for 50 typical cycles (including sterilization and use). The actual number of uses may vary depending on usage conditions. Discard if any End of Life Indicators are present.
- This product is supplied non-sterile and must be sterilized prior to use.

### Indications/Intended Use

The Reusable Articulating Inserts are designed for use with Encision's AEM® Handle Assembly or e-Access™ Slide Lock Handle.

Scissors Inserts are intended for use on *soft* tissue only.

### Contraindications

These inserts are not intended for use when laparoscopic techniques are contraindicated.

### Instructions For Use

#### Prior to Use

Thoroughly read these instructions and the instructions in the Handle Instructions for Use.

Reusable components are supplied non-sterile. Refer to Reprocessing Instructions prior to use.

Inspect the insert and handle for proper assembly and function. See End of Life Indicators.

### CAUTION

- Wear of insert housing may be an indication of impending mechanical failure, which may result in component coming loose during use.

### Assembly Instructions

See Handle Instructions for Use for assembly of insert into handle.

### NOTE

- Used instruments are considered medical waste. Dispose of in accordance with local regulations.

### Reprocessing Instructions

After use, the insert must be cleaned, lubricated and sterilized before reuse.

This product has been validated for 10<sup>-6</sup> sterility assurance level when reprocessed in accordance with these cleaning and sterilization instructions.

### CAUTION

- Treat a used insert as a potential biohazard until cleaning and sterilization has been completed. Microscopic residues may remain after cleaning.
- Stainless Steel is rust-resistant, not rust-proof. Inserts should be kept dry when not in use.
- Do not use bleach (sodium hypochlorite) based products during cleaning. Bleach is extremely corrosive to metals and can negatively affect the electrosurgical instruments. Refer to your cleaning products Material Safety Data Sheet (MSDS) to ensure that they are not corrosive or harmful to various metals (including stainless steel, gold, etc).
- Certain cleaning chemicals may negatively affect metals. Refer to your cleaning products Material Safety Data Sheets (MSDS) to ensure that they are not corrosive or harmful to various metals.

### Cleaning


Any of the three following cleaning methods may be used: Manual Cleaning, Combination Manual/Automated Cleaning, or Automated Cleaning.

### NOTE

- The use of fully distilled or demineralized water is preferred for cleaning and rinsing.

#### Option 1: Manual Cleaning

Step	Description
1	Immerse the insert in KLENZYME® or equivalent blood dissolving enzymatic solution (prepared per manufacturer's instructions) for at least 5 minutes and gently agitate. (Soak longer if proteinaceous material is present.)
2	Remove the insert from the enzyme solution and rinse thoroughly under running tap water (minimum of 3 minutes).
3	Immerse the insert in MANU-KLENZ® or equivalent detergent solution (prepared per manufacturer's instructions) and clean all surfaces.

4	Using a hand held soft bristle brush, use a back and forth motion and brush all surfaces with particular attention paid to crevices, grooves, fittings, and joints. 
5	Remove the insert from the detergent solution and rinse thoroughly under running distilled or demineralized water (minimum of 3 minutes).

#### Option 2: Combination Manual/Automated Cleaning

Step	Description
1	Follow Manual Cleaning steps 1 – 5.
2	Use of a sonicator or ultrasonic cleaner at 35-45 kHz can support manual cleaning of devices. Cycle time 5 minutes maximum. Do not exceed water temperatures above 50°C (122°F).
3	Place inserts in a suitable washer/disinfector basket and process through a standard instrument washer/disinfector cleaning cycle as described below.

#### Option 3: Automated Cleaning

Preparation: Remove coarse debris from the instruments immediately after each use, paying particular attention to the cord connector, crevices, grooves, fittings, joints, and lumens. A soft bristle brush is recommended.

Instruments may be run through a sonicator or ultrasonic cleaner at 35-45 kHz prior to the automated washing cycle listed below. Cycle time 5 minutes maximum. Do not exceed water temperatures above 50°C (122°F)

Step	Description
1	Place the product on a tray that is suitable for cleaning.  If the instrument has a lumen or channel connect it to the water injectors of the cleaning/disinfection machine.  Store any products with hinges or joints on the tray in such a way that the joints are open.
2	Prerinse for a minimum of 3 minutes with cold water.
3	Discharge.
4	Wash for a minimum of 5 minutes at 60± 5°C (140±9°F) with a 0.5% alkaline cleaning agent or an alkaline cleaning agent with a chemistry of 10.9pH.  Water temperatures up to 95°C (203°F) may be used, as recommended by the chemical manufacturer.
5	Discharge.

6	Neutralize as advised by the chemical manufacturer for a minimum of 3 minutes with hot water (>40°C, >104°F) and a neutralizing agent.  Water temperatures up to 95°C (203°F) may be used, as recommended by the chemical manufacturer.
7	Discharge.
8	Rinse for a minimum of 3 minutes with hot water (>40°C, >104°F).  Water temperatures up to 95°C (203°F) may be used.
9	Discharge.
10	Thermal disinfection and drying step for a minimum of 5 minutes at 90°C (194°F). Lubricant may be added during this cycle.  Temperature may be as high as 135°C (275°F).
Drying: If necessary, manual drying may be carried out using a lint free cloth. Dry cavities by blowing with sterile compressed air.	

**NOTE**

- The washer/disinfector manufacturer's instructions should be strictly adhered to.

**Lubrication**

After cleaning all the components or as part of the automated cleaning cycle, lubricate the shaft and all moving parts of the jaw insert with water-soluble medical instrument lubrication. After lubrication, prepare the instrument for sterilization by wrapping in sterile wrap, as appropriate.

**NOTE**

- Failure to lubricate may result in increased wear or binding.

**Sterilization**

Monitoring sterility using Geobacillus stearothermophilus spore strips is recommended.

**NOTE**

- Do not exceed temperatures of 138°C / 280°F. Performance to specification has not been verified above this temperature. Damage to the product may occur.
- Validation of sterilization cycles other than steam autoclave is the responsibility of the end user.

Step	Description
1	Allow the insert to dry thoroughly prior to sterilization.
2	See Caution above regarding removal of detergent solution from the connector surface using an alcohol swab prior to sterilization.
3	Steam Sterilization (use one of the methods listed below):  Prevac (2 layers of Sterile Wrap or Unwrapped)  Temperature: 132°C / 270°F minimum Duration: 4 minutes minimum

**Gravity Steam (2 layers of Sterile Wrap)**

Temperature: 132°C / 270°F minimum  
Duration: 15 minutes minimum.

**Gravity Steam (FLASH/Unwrapped)\***

Temperature: 132°C / 270°F minimum  
Duration: 10 minutes minimum

\* Flash Sterilization is not recommended, except in emergency situations.

**CAUTION**

- Cold Soak Sterilization is not recommended. Testing results indicate that cold soak sterilization (Glutaraldehyde) is not adequate for this product. Manufacturer's recommended cycle times DO NOT always provide a 10<sup>-6</sup> sterility assurance level.

**Product Life**

This product has been tested and shown to remain safe and effective for up to 50 uses under typical use and reprocessing conditions.

The life of this surgical instrument is largely dependent on the care and handling at the point of use / cleaning / sterilization. For optimal instrument life, protect it from contact / impact with other instruments during decontamination and sterilization.

The number of uses may be reduced by improper handling, and "prevac" or "flash" sterilization methods.

Electrification of scissors inserts may reduce their useful life.

**End of Life Indicators**

The inserts may be used until an insert exhibits any of the following end of life indicators.

Visually inspect prior to use. Discontinue use if any of the following are evident:

- Bent shaft or tip affecting insertion into lumen of handle
- Dulling of Scissors
- Binding or impaired mechanical functions
- Wear of housing as evidenced by abrasion of Encision name

**Express Warranty**

Encision hereby warrants to Buyer that products purchased hereunder shall be free from defects in materials and workmanship under normal use and service, as specified in these Instructions for Use/Care, for the period of six months from date of purchase, or upon reaching an end of life indicator (including Maximum Number of Uses), whichever occurs first. Any evidence of unauthorized modification or repair of this device will void the warranty.

**Return of Used Product**




If for any reason this product must be returned to Encision, a returned goods authorization is required prior to shipping. Appropriate return instructions may be obtained from Encision.

**Product**

Encision reserves the right to amend, modify or to change any product, to introduce new products, to withdraw products and otherwise vary product specifications at any time without notice.

Encision®, AEM®, and e-Access™ are trademarks or registered trademarks of Encision Inc.

KLENZYME® and MANU-KLENZ® are registered trademarks of Steris Corporation.

SYMBOL	MEANING	ISO 7000 REF NO.
	Manufacturer	3082
	Date of Manufacture	2497
	Batch Code or Lot Code	2492
	Catalog Number/Reference	2493
	Medical Device that has not been subjected to a sterilization process when similar or identical products are sold Sterile.	2609
	Consult Instructions for Use	1641
	CAUTION - Read Instructions for use for further details	0434A
SYMBOL	MEANING	REFERENCE
	Authorized Representative in the European Community	ISO 15223 5.1.2 EN980 5.13
	CE Mark	Medical Device Directive 93/42/EEC
	Does not contain or no presence of natural rubber latex	ISO 15223 5.4.5 with negation symbol Annex B.2
	Prescription Only	21 CFR PART 801

**Manufactured by**

Encision Inc.  
6797 Winchester Circle  
Boulder, CO 80301 USA  
Ph: 800-998-0986 Fax: 303-444-2693  
www.encision.com info@encision.com

**Made In USA****EC REP Authorized Representative**

(according to MDD93/42/EEC)  
MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany

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© Copyright 2016 Encision Inc.
Encision®, AEM® en e-Access™ zijn handelsmerken of gedeponeerde handelsmerken van Encision Inc.
KLENZYME® en MANU-KLENZ® zijn gedeponeerde handelsmerken van Steris Corporation.